

ctDNA analysis of lymph detects minimal residual disease and correlates with recurrence in HPV-negative head and neck cancer patients



Wendy Winckler¹, Zhuosheng Gu¹, Seka Lazare¹, Adam Benson¹, Noah Earland², Marra S. Francis¹, Adam Harmon¹, Megan Long¹, Maciej Pacula¹, Ashley Tellis¹, Aadel A. Chaudhuri³, Jose P. Zevallos⁴


¹Droplet Biosciences, Cambridge, MA ²Washington University, St. Louis, MO

³Mayo Clinic, Rochester, MN ⁴University of Pittsburgh Medical Center, Pittsburgh, PA


Introduction

Relapse is a major cause of failure after surgery in human papillomavirus (HPV) negative head and neck squamous cell carcinoma (HNSCC), with up to 50% of patients recurring within 2 years. Clinicopathologic criteria for adjuvant treatment are imprecise and have not changed for decades. Using an ultra-sensitive targeted sequencing approach, we demonstrate that circulating tumor DNA (ctDNA) in lymphatic exudate collected via surgical drains ("lymph") 24 hours after surgery identifies MRD and outperforms plasma in an initial cohort of 36 HNSCC patients. We then replicated the lymph ctDNA test in an independent, multi-site cohort of 37 HNSCC patients. Lymph performance was particularly enhanced in locoregional relapse in both cohorts and generalized to early-stage patients. Analysis of matched plasma collected at this early timepoint was not predictive of recurrence. We demonstrate a liquid biopsy approach using a historically overlooked biofluid to potentially enable precision adjuvant therapy and achieve superior oncologic outcomes. (Abstract updated with additional data and analyses.)


Methods and Materials




Lymph and blood were collected from 76 HPV-negative HNSCC patients postoperatively at 24 hours along with resected tumor and sequenced using a custom HNSCC-focused 6.3 Mb panel. Somatic mutations were identified by sequencing tumor and blood (>250x). Lymph and plasma were sequenced to >6,000x coverage. Patient-specific somatic variants detected in tumor were directly genotyped to examine those positions in lymph and plasma using a custom pipeline. One patient was excluded due to a blood sample mismatch and two patients were excluded because tumor samples failed sequencing, resulting in a total study population of 73 patients. Of the 73 patients, 36 patients were studied in an initial cohort and 37 were reserved for the replication cohort. Forty-two patients had matched plasma samples available (29 initial cohort, 13 replication cohort). A sample was classified as positive if the mean ctDNA variant allele fraction (VAF) exceeded a calibrated cut-off accounting for tumor mutation burden and coverage. The Kaplan-Meier estimator with log-rank test and Cox proportional-hazards model were used for survival analyses. Logistic regression models were performed with 5-fold cross-validation.



Swimmer plot of recurred (n=31) and non-recurred (n=42) patients describing lymph ctDNA status, treatment course, recurrence type and outcome.




Lymph but not plasma ctDNA correlates with recurrence




Results



Locoregional recurrences only



ENE-negative patients



Pathology feature performance in this cohort



Conclusions

- Post-surgical lymphatic fluid is a practical, differentiated proximal biofluid for ctDNA MRD detection 24 hours after head and neck cancer surgery
- Lymph ctDNA stratified survival outcomes at the 24h post-surgical timepoint while plasma ctDNA did not, and demonstrated enhanced performance in locoregional recurrence
- Lymph ctDNA significantly stratifies patients who lack high risk pathology features such as extranodal extension, outperforming individual pathology features and synergizing in a combination model
- Detection of lymph ctDNA 24 hours after head and neck cancer surgery could represent a precision biomarker for adjuvant decision-making
- Future studies will validate these findings in larger prospective multi-institutional settings and extend this approach to additional cancer types

Contact Information

- For more information, please contact: Wendy Winckler at wwinckler@dropletbiosci.com